

661 lbs
standard weight capacity

Large battery
for extended performance

STILLE imagiQ2™ surgical imaging table

Setting new low-dose standards in the OR

The STILLE imagiQ2™ surgical imaging table represents an important advance in developing the low-dose technology field, making the OR an even safer place, for patients and for the entire surgical team.

More efficient and shorter procedure times

The True Free Float™ technology of the STILLE imagiQ2 gives unprecedented control and freedom of movement in any desired direction while providing state-of-the-art smoothness. Like an extension of the surgeon's hand, the True Free Float™ increases efficiency in the OR and shortens procedure times.

Choose a Low-Dose Enabler™

The exceptional high-performance carbon fiber table top, provides outstanding radiolucency of 0.4mm Al equivalency, enabling the C-arm to reduce radiation dose without compromising image quality. That is 60% more radiolucent than average tables on the market.

Optimizing C-arm capabilities

With its ultra-slim design and super-low base, the imagiQ2 enhances C-arm access, making it the ideal C-arm complement

for the mobile Hybrid OR or office suite. Versatility through multiple surgical specialties including: Vascular/Cardiac, ERCP/Endoscopy, Pain Management/Neurosurgery and Orthopedics. These features minimize the need for C-arm adjustments that could contribute to overheating.

Continuous image focus

The patented technology of the ISO-ROLL™ cradling motion ($\pm 15^\circ$) helps maintain image focus even while in full lateral roll, minimizing patient and C-arm repositioning, as well as reducing scatter radiation.


Gold-standard service and warranty

STILLE is dedicated to providing gold-standard customer service. We believe in our products: imagiQ2 is backed with a 2-year warranty on labor and 3 years on parts. Some restrictions apply. For full warranty terms please see our website: www.stille.se/us

Table movement ranges⁽¹⁾

Height adjustment	28.0"–43.0"
Trendelenburg	±20°
Lateral roll (ISO-ROLL™)	±15°
Longitudinal True Free Float™	35.4"
Lateral True Free Float™	9.9"
Diagonal True Free Float™	36.8"

⁽¹⁾ ± 2.0° tolerance for angles, ± 0.5" for lengths



Product specifications of the STILLE imagiQ2™

Part number	530-1703
Weight capacity	661 lbs in any position; the table top does not have to be centered in lock mode.
Table top length, width	L= 90.5" ; W= 21.6"
Radiolucent length	89.0" with 3D headrest / 66.6" without headrest
Float & lock	True Free Float™ allowing smooth movements + STILLE quick-lock system
AC power: 100–250 Volts	Direct AC power + Power Performance Battery DC power
Side rails	Fixed, about 56.1"
Exchangeable table top	Yes, future-proof platform for extended use
Standard configuration	Includes: 2 hand controls, 1 extension cable, 1 pan handle, mattress

Available accessories⁽²⁾ ⁽³⁾

Top individual best-sellers

- ▶ 3D head rest (535-1745)
- ▶ integrated control module US (535-1766)
- ▶ hand control holder (514-50-3SR040)
- ▶ adjustable head rest (535-1742)
- ▶ foot control (535-1762)
- ▶ flexible anesthesia screen (512-5030)
- ▶ IV pole without clamp (512-5010)
- ▶ IV pole clamp (512-5013)

Accessory package⁽⁴⁾ (part number: 530-1790)

- ▶ long detachable side rail (535-1720)
- ▶ short detachable side rail (535-1724)
- ▶ fistula arm board (535-1730)
- ▶ catheter tray 39.4" (535-1750)
- ▶ radiation shield loop, grey (535-1772)
- ▶ deluxe patient strap (O-RSX2)

⁽²⁾ 1-year warranty on all accessories.

⁽³⁾ See the accessories brochure for the full list.

⁽⁴⁾ Package offers discount over individual items purchased separately. Items of this package can be ordered individually.

STILLE AB
Ekbacken 11
644 30 Torshälla
Sweden

☎ +46 8 588 580 00
✉ sales@stille.se
🌐 www.stille.se

STILLE Surgical Inc.
840 Parkview Blvd
Lombard, IL 60148
United States

☎ 1.800.665.1614
✉ sales.us@stille.se
🌐 www.stille.se/us

These products may not be available in all regions. Please contact your STILLE-representative for more information on market availability.

This document correctly describes these devices at the time the document was printed. Please always refer to the user manual for any product amendments.

The product images shown are for illustration purposes only and may not be an exact representation of the product.